
TRPM6 Antibody Blocking Peptide

Catalog Number: bs-9048P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: TRPM6 is an essential ion channel and serine/threonine-protein kinase, and is crucial for magnesium homeostasis. TRPM6 also has an important role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. The various isoforms of the type M6-kinase lack the ion channel region.