

---

## TENM3 Antibody Blocking Peptide

Catalog Number: bs-9061P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of the teneurin transmembrane protein family. The encoded protein may be involved in the regulation of neuronal development including development of the visual pathway. Mutations in this gene have been associated with microphthalmia and developmental dysplasia of the hip. [provided by RefSeq, Jan 2023]