

---

## apoptin Rabbit pAb

Catalog Number: bs-8720R

Target Protein: apoptin

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Chicken anemia virus)

Predicted MW: 13 kDa

Source: KLH conjugated synthetic peptide derived from Chicken anemia virus apoptin: 51-121/121.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Apoptin, a protein encoded by vp3 gene of chicken anemia virus, can specifically induce apoptosis of tumor cells but has no toxic effect on normal cells.