

bs-8720R

[Primary Antibody]

apoptin Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>Target: apoptin</p> <p>Immunogen: KLH conjugated synthetic peptide derived from Chicken anemia virus apoptin: 51-121/121.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Apoptin, a protein encoded by vp3 gene of chicken anemia virus, can specifically induce apoptosis of tumor cells but has no toxic effect on normal cells.</p>	<p>Isotype: IgG</p> <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Chicken anemia virus)</p> <p>Predicted MW.: 13 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	--