

---

## TRIM35 Antibody Blocking Peptide

Catalog Number: bs-9149P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: TRIM35 is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1, a B-box type 2 and a coiled-coil region. TRIM35 may play a role as a tumor suppressor and is implicated in the cell death mechanism. There are two named isoforms.