

bs-8660R

[Primary Antibody]

Silent protein UshA(0) Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: SALTY) Predicted MW.: 58 kDa
Clonality: Polyclonal		
Target: Silent protein UshA(0)		
Immunogen: KLH conjugated synthetic peptide derived from SALTY of Salmonella typhimurium (O): 301-400/550.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=26.8]** Chen Yanmei. et al. Reinforcement of the intestinal mucosal barrier via mucus-penetrating PEGylated bacteria. NAT BIOMED ENG. 2024 Jun;;1-19 IF ;Mouse. 38839928
- **[IF=26.8]** Guocheng Xie. et al.Enhancing phage therapy by coating single bacteriophage-infected bacteria with polymer to preserve phage vitality.NATURE BIOMEDICAL ENGINEERING.2025 Feb 25. IF ;Mouse. 40000896