
RNF34 Antibody Blocking Peptide

Catalog Number: bs-9242P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Has E3 ubiquitin-protein ligase activity. Regulates the levels of CASP8 and CASP10 by targeting them for proteasomal degradation. Protects cells against apoptosis induced by TNF. Binds phosphatidylinositol-5-phosphate and phosphatidylinositol-3-phosphate.