

bs-19022R

[Primary Antibody]

NARG1L Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 79612 Target: NARG1L Immunogen: KLH conjugated synthetic peptide derived from human NARG1L: 551-650/864. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: The NARG1-like protein (NARG1L) is an 864 amino acid protein that contains 7 TPR repeats and may also be a component of a complex that displays N-terminal acetyltransferase activity.	Isotype: IgG SWISS: B0QZ59	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Cow, Chicken, Dog, Horse) Predicted MW.: 101 kDa Subcellular Location: Cytoplasm
--	--	--