

---

## **BCL7B Antibody Blocking Peptide**

Catalog Number: bs-9298P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: BCL7B shows high homology to the BCL7A protein, which is known to be directly involved in a complex chromosomal translocation in Burkitt lymphoma cell lines. The specific function of BCL7B has not yet been determined. However, it may play a role in lung tumor development or progression. The BCL7B gene is located at a chromosomal region commonly deleted in the congenital disorder, Williams syndrome.