
RRP1B Antibody Blocking Peptide

Catalog Number: bs-9300P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: RRP15 is a 282 amino acid protein that localizes to the nucleus and may be involved in the modification and processing of precursor rRNA, possibly playing a role in maturation of the 60S ribosomal subunit. The gene encoding RRP15 maps to human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.