
PGBD1 Antibody Blocking Peptide

Catalog Number: bs-9370P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: PGDB1, also known as SCAND4 or HUCEP-4, is an 809 amino acid protein that contains one SCAN box domain and is a member of the PGBD family. Expressed specifically in brain tissue, PGBD1 may, via its SCAN box domain, be involved in transcriptional regulation events within the nucleus.