
MYCBP2 Antibody Blocking Peptide

Catalog Number: bs-9375P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: MYCBP2 belongs to the highwire family. It is a probable E3 ubiquitin-protein ligase which mediates ubiquitination and subsequent proteasomal degradation of target proteins. MYCBP2 may function as a facilitator or regulator of transcriptional activation by MYC and have a role during synaptogenesis. There are two different isoforms.