
TSPAN17 Antibody Blocking Peptide

Catalog Number: bs-9424P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: TSPAN17 belongs to the tetraspanin (TM4SF) family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. There are four named isoforms.