

---

## ALDOB Antibody Blocking Peptide

Catalog Number: bs-4026P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ALDOB is a tetrameric glycolytic enzyme that catalyzes the reversible conversion of fructose-1,6-bisphosphate to glyceraldehyde 3-phosphate and dihydroxyacetone phosphate. Defects in ALDOB cause hereditary fructose intolerance.