

## Recombinant human TNFSF13B protein, His (HEK293)

Catalog Number: bs-43027P

Concentration: >0.5 mg/ml

AA Seq: 134-285/285

Predicted MW: 21

Detected MW: 21 kDa

Tags: N-His

Activity: Yes

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B/TACI, TNFRSF17/BCMA, and TNFRSF13C/BAFFR. This cytokine is expressed in B cell lineage cells, and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. [provided by RefSeq].

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.