
Human IgG FC fragment, Biotin conjugated

Catalog Number: bs-0365P-Bio

Concentration: 1.0 mg/ml

AA Seq: Purified native protein

Activity: Not tested

Endotoxin: Not analyzed

Form: Liquid

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.