
hnRNP F Antibody Blocking Peptide

Catalog Number: bs-4205P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: hnRNP F is part of the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). These are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA) and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. hnRNP F binds to RNAs which have guanosine-rich sequences.