

bs-6501R

[Primary Antibody]

Rab25 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 57111</p> <p>Target: Rab25</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human Rab25: 145-190/213.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: RAB proteins, such as RAB25, are members of the RASsuperfamily of small GTPases that are involved in membranetrafficking. Members of the RAB11 subfamily, including RAB25,control the return of internalized membrane-associated moieties tothe cell surface</p>	<p>Isotype: IgG</p> <p>SWISS: P57735</p>	<p>Applications: ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Horse)</p> <p>Predicted MW.: 23 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
--	--	---