
Bovine Serum Albumin, BF647 conjugated

Catalog Number: bs-0292P-BF647

Concentration: 1mg/ml

AA Seq: Purified native protein

Activity: Not tested

Endotoxin: Not analyzed

Form: Liquid

Storage: 10 mM TBS (pH=7.4) with 0.02% Proclin300 and 50% glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: Albumin is one of the most extensively utilized proteins in the medical industry. Its availability, inherent stability, and ligand binding ability make albumin an essential component in the production of many therapeutic drugs, vaccines and diagnostic tests.