
CCAP Antibody Blocking Peptide

Catalog Number: bs-9859P

Activity: Yes

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Cardioactive peptide (CCAP) is a cardioregulatory neurohormone that increases heart beat rate during adult wing inflation. This protein has no effect on beat amplitude. The effect of cardioactive peptide is both ino- and conotrophic. Cardioactive peptide has tissue specificity within the central nervous system and most neurons exhibit coexpression with burs.