

---

## **CRHBP Antibody Blocking Peptide**

Catalog Number: bs-9906P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: CRHBP is a potent stimulator of synthesis and secretion of preopiomelanocortin derived peptides. Although CRH concentrations in the human peripheral circulation are normally low, they increase throughout pregnancy and fall rapidly after parturition. Maternal plasma CRH probably originates from the placenta. Human plasma contains a CRH binding protein which inactivates CRH and which may prevent inappropriate pituitary adrenal stimulation in pregnancy.