

---

## **SPINT1 Antibody Blocking Peptide**

Catalog Number: bs-9913P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: HGF activator belongs to the peptidase family S1 and activates hepatocyte growth factor (HGF) by converting it from a single chain to a heterodimeric form. HGF Activator is produced as a zymogen, and the form found in serum is largely the inactive single-chain form. It is activated to a heterodimeric form by endoproteolytic processing.