

---

## Marijuana/BSA

Catalog Number: bs-2028PB

Concentration: 1mg/ml

Activity: No

Form: Lyophilized or Liquid

Storage: 0.01M TBS(pH7.4).

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Cannabinol (CBN) is a weakly psychoactive cannabinoid found in Cannabis sativa and Cannabis indica. It is also a metabolite of tetrahydrocannabinol (THC). CBN acts as a weak agonist of CB1 receptors but has a higher affinity to CB2 receptors, with lower affinities in comparison to THC. Because of its somewhat selective CB2 agonism, it is used experimentally as an immunosuppressant.