bs-10630R

[Primary Antibody]

Bioss ANTIBODIES

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

– DATASHEET –

Host: Rabbit Isotype: IgG

Clonality: Polyclonal
GenelD: 100502522
Target: ZmCCD8

ZmCCD8 Rabbit pAb

Immunogen: KLH conjugated synthetic peptide derived from ZmCCD8:

441-540/572.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: Zea mays)

Predicted MW.: 63 kDa