
MRAP2 Rabbit pAb

Catalog Number: bs-11419R

Target Protein: MRAP2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 24 kDa

Subcellular: Cell membrane ,Cytoplasm

Locations:

Entrez Gene: 112609

Swiss Prot: Q96G30

Source: KLH conjugated synthetic peptide derived from human MRAP2: 81-130/205.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: No data available.