

bs-11419R

[Primary Antibody]

MRAP2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 112609	SWISS: Q96G30	Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)
Target: MRAP2		
Immunogen: KLH conjugated synthetic peptide derived from human MRAP2: 81-130/205.		
Purification: affinity purified by Protein A		Predicted MW.: 24 kDa
Concentration: 1mg/ml		Subcellular Location: Cell membrane ,Cytoplasm
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: No data available.		