
Recombinant EBFP protein, His

Catalog Number: bs-33185P

Concentration: >0.5 mg/ml

AA Seq: 1-238/238

Predicted MW: 26

Detected MW: 26 kDa

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Storage: Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.