

bs-4477R

[Primary Antibody]

PLEKHA9 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 51054 Target: PLEKHA9 Immunogen: KLH conjugated synthetic peptide derived from human PLEKHA9: 101-200/391. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: PLEKHA8P1 (Pleckstrin Homology Domain Containing A8 Pseudogene 1) is a Pseudogene. GO annotations related to this gene include glycolipid binding and glycolipid transporter activity.	Isotype: IgG SWISS: O95397	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human) Predicted MW.: 43 kDa Subcellular Location: Cytoplasm
--	--	---