

bs-12728R

[Primary Antibody]

PLET1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500)
Clonality: Polyclonal		IHC-F (1:100-500)
GeneID: 349633	SWISS: Q6UQ28	IF (1:100-500)
Target: PLET1		ICC/IF (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from human PLET1: 41-140/204.		Reactivity: (predicted: Human)
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		Predicted MW.: 23 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Subcellular Location: Cell membrane