

Human FSH Mouse mAb

Cat. No: V6505
Name: Human FSH Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5H7
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Concentration: >0.5mg/ml
Buffer: 0.01M PBS (pH7.4).
Antigen: Recombinant human FSH protein ([bs-43084P](#))
Pairing Information:

ELISA		Detection				
		V6501	V6502	V6503	V6504	V6505
Capture	V6501	-	-	-	-	-
	V6502	-	-	-	-	-
	V6503	-	-	-	-	+++
	V6504	-	-	-	-	+++
	V6505	-	-	-	-	-
ICA		Detection				
		V6501	V6502	V6503	V6504	V6505
Capture	V6501	-	-	-	-	+++
	V6502	-	-	-	-	+
	V6503	-	-	-	-	+
	V6504	-	-	-	-	-
	V6505	-	-	-	-	-