

Human FSH Mouse mAb

Cat. No: V6503
Name: Human FSH Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 28H1
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Lyophilized or Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human FSH Protein ([bs-43084P](#))
Pairing Information:

ELISA		Detection				
		V6501	V6502	V6503	V6504	V6505
Capture	V6501	-	-	-	-	-
	V6502	-	-	-	-	-
	V6503	-	-	-	-	+++
	V6504	-	-	-	-	+++
	V6505	-	-	-	-	-
ICA		Detection				
		V6501	V6502	V6503	V6504	V6505
Capture	V6501	-	-	-	-	+++
	V6502	-	-	-	-	+
	V6503	-	-	-	-	+
	V6504	-	-	-	-	-
	V6505	-	-	-	-	-