

Human HCG-α Mouse mAb

Cat. No: V2702
Name: Human HCG-α Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 2A5
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ICA
Antigen: Recombinant Human HCG-α Protein ([bs-43085P](#))
Pairing Information:

ELISA		Detection						
		V2701	V2702	V2703	V2704	V2705	V2706	V2707
Capture	V2701	-	-	-	-	-	-	-
	V2702	-	-	-	-	-	-	-
	V2703	-	-	-	-	-	-	-
	V2704	-	-	-	-	-	-	-
	V2705	-	-	-	-	-	-	+++
	V2706	-	-	-	-	-	-	+++
	V2707	-	-	-	-	-	-	-
ICA		Detection						
		V2701	V2702	V2703	V2704	V2705	V2706	V2707
Capture	V2701	-	-	-	-	-	-	+
	V2702	-	-	-	-	-	-	+++
	V2703	-	-	-	-	-	-	+
	V2704	-	-	-	-	-	-	+
	V2705	-	-	-	-	-	-	-
	V2706	-	-	-	-	-	-	-
	V2707	-	-	-	-	-	-	-