

## Human HCG- $\alpha$ Mouse mAb

|                      |                                                                                 |
|----------------------|---------------------------------------------------------------------------------|
| Cat. No:             | V2701                                                                           |
| Name:                | Human HCG- $\alpha$ Mouse mAb                                                   |
| Storage:             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:          | Monoclonal                                                                      |
| Clone No:            | 7D11                                                                            |
| Host/Source:         | Mouse                                                                           |
| Subtype:             | IgG1                                                                            |
| Reactivity:          | (predicted:&nbsp;Human)                                                         |
| Format:              | Liquid                                                                          |
| Buffer:              | 0.01M PBS (pH7.4).                                                              |
| Applications:        | ICA                                                                             |
| Antigen:             | Recombinant Human Human HCG- $\alpha$ Protein ( <a href="#">bs-43085P</a> )     |
| Pairing Information: |                                                                                 |

| ELISA   |       | Detection |       |       |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|-------|-------|
|         |       | V2701     | V2702 | V2703 | V2704 | V2705 | V2706 | V2707 |
| Capture | V2701 | -         | -     | -     | -     | -     | -     | -     |
|         | V2702 | -         | -     | -     | -     | -     | -     | -     |
|         | V2703 | -         | -     | -     | -     | -     | -     | -     |
|         | V2704 | -         | -     | -     | -     | -     | -     | -     |
|         | V2705 | -         | -     | -     | -     | -     | -     | +++   |
|         | V2706 | -         | -     | -     | -     | -     | -     | +++   |
|         | V2707 | -         | -     | -     | -     | -     | -     | -     |
| ICA     |       | Detection |       |       |       |       |       |       |
|         |       | V2701     | V2702 | V2703 | V2704 | V2705 | V2706 | V2707 |
| Capture | V2701 | -         | -     | -     | -     | -     | -     | +     |
|         | V2702 | -         | -     | -     | -     | -     | -     | +++   |
|         | V2703 | -         | -     | -     | -     | -     | -     | +     |
|         | V2704 | -         | -     | -     | -     | -     | -     | +     |
|         | V2705 | -         | -     | -     | -     | -     | -     | -     |
|         | V2706 | -         | -     | -     | -     | -     | -     | -     |
|         | V2707 | -         | -     | -     | -     | -     | -     | -     |