

Human HCG-beta Mouse mAb

Cat. No:	V2707
Name:	Human HCG-beta Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4H7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant HCG-beta Protein (bs-43172P)
Pairing Information:	

ELISA		Detection						
		V2701	V2702	V2703	V2704	V2705	V2706	V2707
Capture	V2701	-	-	-	-	-	-	-
	V2702	-	-	-	-	-	-	-
	V2703	-	-	-	-	-	-	-
	V2704	-	-	-	-	-	-	-
	V2705	-	-	-	-	-	-	+++
	V2706	-	-	-	-	-	-	+++
	V2707	-	-	-	-	-	-	-
ICA		Detection						
		V2701	V2702	V2703	V2704	V2705	V2706	V2707
Capture	V2701	-	-	-	-	-	-	+
	V2702	-	-	-	-	-	-	+++
	V2703	-	-	-	-	-	-	+
	V2704	-	-	-	-	-	-	+
	V2705	-	-	-	-	-	-	-
	V2706	-	-	-	-	-	-	-
	V2707	-	-	-	-	-	-	-