

Recombinant human IL-10 protein, N-His-Avi (HEK293), Biotin conjugated

Catalog Number: bs-47119P-Bio

Concentration: >0.5 mg/ml

AA Seq: 19-178/178

Predicted MW: 17

Detected MW: Due to glycosylation, the protein migrates to 23-25 kDa based on Tris-Bis PAGE result.

Tags: N-His-Avi

Activity: Not tested

Endotoxin: <1.0 EU/μg as determined by LAL

Purity: >95% as determined by Tris-Bis PAGE


Purification: AC

Storage: Lyophilized from 0.22um filtered solution in PBS (pH7.4) with 5mM DTT. Normally 5% trehalose is added as protectant before Lyophilization.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: The protein encoded by this gene is a cytokine produced primarily by monocytes and to a lesser extent by lymphocytes. This cytokine has pleiotropic effects in immunoregulation and inflammation. It down-regulates the expression of Th1 cytokines, MHC class II Ags, and costimulatory molecules on macrophages. It also enhances B cell survival, proliferation, and antibody production. This cytokine can block NF-kappa B activity, and is involved in the regulation of the JAK-STAT signaling pathway. Knockout studies in mice suggested the function of this cytokine as an essential immunoregulator in the intestinal tract. Mutations in this gene are associated with an increased susceptibility to HIV-1 infection and rheumatoid arthritis. [provided by RefSeq, May 2020]

VALIDATION IMAGES



Human IL-10 Protein on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.