
Orai1 Antibody Blocking Peptide

Catalog Number: bs-18019P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Calcium release activated calcium channel protein 1; Calcium release activated calcium modulator 1; Calcium release-activated calcium channel protein 1; CRACM 1; CRACM1; CRACM1; CRACM1_HUMAN; FLJ14466; Orai 1; ORAI calcium release activated calcium modulator 1; orai1; ORAT1; Protein orai 1; Protein orai-1; TMEM 142A; TMEM142A; Transmembrane protein 142A.