

Recombinant human Trem2 protein, C-His-Avi (HEK293), Biotin conjugated

Catalog Number: bs-47221P-Bio

Concentration: >0.5 mg/ml

AA Seq: 19-174/230

Predicted MW: 17

Detected MW: Due to glycosylation, the protein migrates to 30-45 kDa based on Tris-Bis PAGE result.

Tags: C-His-Avi

Activity: Not tested

Endotoxin: <1.0 EU/μg as determined by LAL

Purity: >95% as determined by Tris-Bis PAGE; >95% as determined by SEC-HPLC


Purification: AC

Storage: Lyophilized from 0.22um filtered solution in PBS (pH7.4) with 5mM DTT. Normally 5% trehalose is added as protectant before Lyophilization.


Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a membrane protein that forms a receptor signaling complex with the TYRO protein tyrosine kinase binding protein. The encoded protein functions in immune response and may be involved in chronic inflammation by triggering the production of constitutive inflammatory cytokines. Defects in this gene are a cause of polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (PLOS). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2012]

VALIDATION IMAGES



Recombinant human Trem2 Protein on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



Recombinant human Trem2 Protein on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.