

**bs-8033R**

**[ Primary Antibody ]**

## COLEC10 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 10584</p> <p><b>Target:</b> COLEC10</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human COLEC10(liver): 151-230/277.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> Lectin that binds to various sugars: galactose &gt; mannose = fucose &gt; N-acetylglucosamine &gt; N-acetylgalactosamine. Tissue specificity:</p>	<p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q9Y6Z7</p>	<p><b>Applications:</b> <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human, Mouse, Rat, Pig, Cow, Dog, Horse)</p> <p><b>Predicted MW.:</b> 28 kDa</p> <p><b>Subcellular Location:</b> Secreted</p>
--	--	---