

bsm-34233A

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rat Anti-Mouse IgD

— DATASHEET —

Host: Rat **Isotype:** IgG2a, κ

Clonality: Monoclonal

Target: Rat Anti-Mouse IgD

Purification: affinity purified by Protein A

Concentration: 0.5mg/ml

Storage: 0.01M PBS (pH7.4).
stored undiluted between 2°C and 8°C.

Applications: Isotype Control
Blocking Assay
etc.
Conjugate-Dependent.

Reactivity: Mouse