

bsm-34233A

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rat Anti-Mouse IgD

— DATASHEET —

Host: Rat	Isotype: IgG2a, κ	Applications: Flow-Cyt (1:100-500)
Clonality: Monoclonal		
Target: Rat Anti-Mouse IgD		
Purification: affinity purified by Protein A		
Concentration: 0.5mg/ml		
Storage: 0.01M PBS (pH7.4). stored undiluted between 2°C and 8°C.		Reactivity: Mouse