

## Recombinant mouse SPRR3 protein, N-His

Catalog Number: bs-42155P

Concentration: 1mg/ml

Species: Mouse

AA Seq: 56-238/238

Predicted MW: 21.5

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).with 150mM NaCl

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** Predicted to enable structural molecule activity. Predicted to be involved in wound healing. Located in Golgi apparatus and perinuclear region of cytoplasm. [provided by Alliance of Genome Resources, Apr 2022]

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.