
Recombinant human AMH protein, C-His

Catalog Number: bs-42006P

Concentration: >0.5 mg/ml

AA Seq: 453-560/560

Predicted MW: 17.2

Detected MW: 20 kDa

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 4M Urea (pH9.0) with 1mM DDT.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Anti-Mullerian hormone is a member of the transforming growth factor-beta gene family which mediates male sexual differentiation. Anti-Mullerian hormone causes the regression of Mullerian ducts which would otherwise differentiate into the uterus and fallopian tubes. Some mutations in the anti-Mullerian hormone result in persistent Mullerian duct syndrome. [provided by RefSeq, Jul 2008].

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.