

---

## Recombinant human ESM1 protein, N-His

Catalog Number: bs-42106P

Concentration: >1mg/ml

AA Seq: 22-184/184

Predicted MW: 19.2

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 2M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** ESM1 is a secreted dermatan sulfate proteoglycan which is predominantly expressed in the endothelial cells in human lung and kidney tissues. Expression of ESM1 is regulated by cytokines. ESM1 may have a role in lung endothelial cell-leukocyte interactions. The transcript contains multiple polyadenylation and mRNA instability signals. Two transcript variants encoding different isoforms have been found for this gene.