
GABA/BSA

Catalog Number: bs-2252PB

AA Seq: Coupling protein

Tags: N/A

Activity: No

Endotoxin: Not analyzed

Form: Liquid

Storage: 0.01M TBS(pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Gamma-aminobutyric acid (GABA) is a major inhibitory neurotransmitter. GABA acts at inhibitory synapses in the brain and spinal cord. Inhibition is provoked by GABA binding resulting in hyperpolarization of the synaptic transmembrane potential of the affected neuron. GABA binding causes ion channels to open allowing either the flow of chloride or potassium ions into or out of the cell.