
Antibacterial protein PR-39 Peptide

Catalog Number: Y-0460

AA Seq: RRRPRPPYLPRRPPPPFFPPRLPPRIPPGFPPRFPPRFP

Predicted MW: 4.72

Purification: HPLC

Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: PR 39 (porcine) is an antibacterial peptide. Stimulates angiogenesis and inhibits inflammatory responses by selectively blocking proteasome degradation of IκBα.