

## Ovalbumin, APC conjugated

Catalog Number: bs-0283P-APC

Concentration: 1mg/ml

AA Seq: Purified native protein

Activity: Not tested


Endotoxin: Not analyzed

Form: Liquid

Storage: 10 mM TBS (pH=7.4) with 0.02% Proclin300 and 50% glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.