

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Ovalbumin, APC conjugated

Catalog Number: bs-0283P-APC

Concentration: 1mg/ml

AA Seq: Purified native protein

Activity: Not tested Endotoxin: Not analyzed


Form: Liquid

Applications: 其它

Storage: 10 mM TBS (pH=7.4) with 0.02% Proclin300 and 50% glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.