

bs-0171R

[Primary Antibody]

DENV-2 E Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: DENV-2 E

Immunogen: KLH conjugated synthetic peptide derived from Dengue virus type-2 envelope glycoprotein E: 501-600/3390.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:500-2000)

ELISA (1:5000-10000)

Reactivity: (predicted: Dengue virus type-2)

Predicted MW.: 54 kDa