

Recombinant mouse FCGR1A (CD64) protein, N-His

Catalog Number: bs-41395P

Concentration: >0.5 mg/ml

AA Seq: 25-297/404

Predicted MW: 30

Detected MW: 33 kDa

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 4M Urea (pH9.0) with 1mM DTT.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1. [provided by RefSeq, Jul 2008]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.