

Recombinant mouse CD68 protein, N-GST & C-His

Catalog Number: bs-42013P

Concentration: >0.5 mg/ml

AA Seq: 119-291/326

Predicted MW: 46.8

Detected MW: 46 kDa

Tags: N-GST & C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 4 M Urea (pH=9.2) with 1 mM DDT.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a 110-kD transmembrane glycoprotein that is highly expressed by human monocytes and tissue macrophages. It is a member of the lysosomal/endosomal-associated membrane glycoprotein (LAMP) family. The protein primarily localizes to lysosomes and endosomes with a smaller fraction circulating to the cell surface. It is a type I integral membrane protein with a heavily glycosylated extracellular domain and binds to tissue- and organ-specific lectins or selectins. The protein is also a member of the scavenger receptor family. Scavenger receptors typically function to clear cellular debris, promote phagocytosis, and mediate the recruitment and activation of macrophages. Alternative splicing results in multiple transcripts encoding different isoforms. [provided by RefSeq, Jul 2008]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.